

OPTETRAK[®]

Operative Technique



OPTETRAK LOGIC[®] ROTATING BEARING KNEE (RBK[®])

.....
Addendum to the Optetrak
Logic LPI Distal First Operative
Technique

INTRODUCTION

Optetrak Logic® is an advanced approach to total knee replacement that includes modern design features and intuitive instrumentation while building on the wisdom of a strong design lineage. With Optetrak Logic RBK®, it's the best of both worlds – the benefits of a fixed and a mobile bearing knee in one unique system.

DESIGN SPECIFICATIONS

The fixed, proximal bearing surface increases flexion potential while maintaining Optetrak's proven optimized congruency. The lower bearing features a patented "wave" design for rotational freedom and predictable kinematics. These two articulating surfaces work together to maintain proper alignment and minimize contact stress.

Follow the Optetrak Logic LPI Distal First Technique for preparation of the femur, tibia and patella. Insert the appropriate Optetrak trials for the final prosthesis trial check.



Figure 1
Place Femoral Trial




Figure 2
Tibial Trial Insert Assembly

Step 1: Place the Optetrak Logic PS femoral trial on the distal femur with the **LPI Locking Femoral Impactor** (*Figure 1*). Ensure that the femoral component is properly positioned on the distal femoral condyles in the medial and lateral direction. Slight upward pressure should be applied to the impactor handle as the component is being impacted to prevent the femoral component from rotating into flexion. Once correct positioning is assured, the component should be fully seated by striking the Locking Femoral Impactor with a mallet.

Place the trial tibial tray and perform trial reduction to assess the stability of the RBK tibial trial insert. If the knee is loose in extension and flexion, proceed to exchange the insert trial with greater thickness and reassess stability.

Step 2: Complete the alignment, stability and patellar tracking checks. Proceed to implantation of final components.



Figure 3
Place RBK Tibial
Prosthesis

Optional: For additional support, screw the **Trial Modular Peg** into the bottom side of the previously selected **RBK Tibial Trial Insert** (*Figure 2*).

Step 3: After the final bone preparation, proceed to the Implantation of Modular Tibial Component, including cementation. Introduce the RBK tibial trays onto the prepared tibial surface and impact using the **Non-Locking RBK Tibial Impactor** (*Figure 3*).

The extraneous cement must be removed from the borders of the tibial component, starting posteriorly and working around to the sides and front. All cement must be removed from the posterior capsular area of the knee. A Tibial Insert Trial should be used when pressurizing the cement during polymerization.

Step 4: Continue the final preparation of the femur and patella as described in the Logic LPI Distal First Technique. Place the RBK polyethylene insert into the previously implanted RBK tibial trial (*Figure 4*). Ensure that the tibial insert is fully seated in the metal tibial tray and proceed to final check and closure.



Figure 4
Implant RBK Tibial Insert

INSTRUMENT LISTING

Catalog Number Part Description

213-65-06 Logic RBK Tibial Impactor



02-013-38-0000 Logic RBK Insert Modular Trial Peg



02-013-38-1009 Logic Tibial Insert Trial, RBK, Size 1, 9mm
 02-013-38-1011 Logic Tibial Insert Trial, RBK, Size 1, 11mm
 02-013-38-1013 Logic Tibial Insert Trial, RBK, Size 1, 13mm
 02-013-38-1015 Logic Tibial Insert Trial, RBK, Size 1, 15mm
 02-013-38-1509 Logic Tibial Insert Trial, RBK, Size 1.5, 9mm
 02-013-38-1511 Logic Tibial Insert Trial, RBK, Size 1.5, 11mm
 02-013-38-1513 Logic Tibial Insert Trial, RBK, Size 1.5, 13mm
 02-013-38-1515 Logic Tibial Insert Trial, RBK, Size 1.5, 15mm
 02-013-38-2009 Logic Tibial Insert Trial, RBK, Size 2, 9mm
 02-013-38-2011 Logic Tibial Insert Trial, RBK, Size 2, 11mm
 02-013-38-2013 Logic Tibial Insert Trial, RBK, Size 2, 13mm
 02-013-38-2015 Logic Tibial Insert Trial, RBK, Size 2, 15mm
 02-013-38-2509 Logic Tibial Insert Trial, RBK, Size 2.5, 9mm
 02-013-38-2511 Logic Tibial Insert Trial, RBK, Size 2.5, 11mm
 02-013-38-2513 Logic Tibial Insert Trial, RBK, Size 2.5, 13mm
 02-013-38-2515 Logic Tibial Insert Trial, RBK, Size 2.5, 15mm
 02-013-38-3009 Logic Tibial Insert Trial, RBK, Size 3, 9mm
 02-013-38-3011 Logic Tibial Insert Trial, RBK, Size 3, 11mm
 02-013-38-3013 Logic Tibial Insert Trial, RBK, Size 3, 13mm
 02-013-38-3015 Logic Tibial Insert Trial, RBK, Size 3, 15mm
 02-013-38-3509 Logic Tibial Insert Trial, RBK, Size 3.5, 9mm
 02-013-38-3511 Logic Tibial Insert Trial, RBK, Size 3.5, 11mm
 02-013-38-3513 Logic Tibial Insert Trial, RBK, Size 3.5, 13mm
 02-013-38-3515 Logic Tibial Insert Trial, RBK, Size 3.5, 15mm
 02-013-38-4009 Logic Tibial Insert Trial, RBK, Size 4, 9mm
 02-013-38-4011 Logic Tibial Insert Trial, RBK, Size 4, 11mm
 02-013-38-4013 Logic Tibial Insert Trial, RBK, Size 4, 13mm
 02-013-38-4015 Logic Tibial Insert Trial, RBK, Size 4, 15mm
 02-013-38-5009 Logic Tibial Insert Trial, RBK, Size 5, 9mm
 02-013-38-5011 Logic Tibial Insert Trial, RBK, Size 5, 11mm
 02-013-38-5013 Logic Tibial Insert Trial, RBK, Size 5, 13mm
 02-013-38-5015 Logic Tibial Insert Trial, RBK, Size 5, 15mm



Exactech is proud to have offices and distributors around the globe.

For more information about Exactech products available in your country, please visit www.exac.com

For additional device information, refer to the Exactech Optetrak Knee System—Instructions for Use for a device description, indications, contraindications, precautions and warnings. • For further product information, please contact Customer Service, Exactech, Inc., 2320 NW 66th Court, Gainesville, Florida 32653-1630, USA. (352) 377-1140, (800) 392-2832 or FAX (352) 378-2617.

This Logic RBK Operative Technique has been developed in cooperation with Michele Brax, MD, Haguenau, France.

Exactech, as the manufacturer of this device, does not practice medicine, and is not responsible for recommending the appropriate surgical technique for use on a particular patient. These guidelines are intended to be solely informational and

each surgeon must evaluate the appropriateness of these guidelines based on his or her personal medical training and experience. Prior to use of this system, the surgeon should refer to the product package insert for comprehensive warnings, precautions, indications for use, contraindications and adverse effects.

The products discussed herein may be available under different trademarks in different countries. All copyrights, and pending and registered trademarks, are property of Exactech, Inc. This material is intended for the sole use and benefit of the Exactech sales force and physicians. It should not be redistributed, duplicated or disclosed without the express written consent of Exactech, Inc. ©2011 Exactech, Inc.

352-377-1140
 1-800-EXACTECH
www.exac.com



712-27-30
 Logic RBK Op.Tech. Addendum 1111

